


CPMOTION™

UNIQUE INNOVATIONS & PRACTICAL BENEFITS

Continuous Passive Motion device is used during the first phase of rehabilitation following lower limb joint surgeries or traumas. The therapy ensures safe exercise by providing passive movement, reduces postoperative pain and minimizes chances of inflammation.



HIGH-SPEED MOVEMENT

- Faster movements for improved therapy outcomes and/or decrease of therapy time



ANKLE MODULE

- Ultimate all-in-one solution – single unit for all 3 leg joints



QUICK AND EASY SET-UP

- Unique click system – no manual tightening needed



REMOTE CONTROL WITH TOUCH SCREEN

- Large display shows 5 parameters on one screen





REAL-TIME ANGLE SETTINGS

- Precise real-time patient's threshold angle setting



ADVANCED THERAPY SETTINGS

- 13 automatic pre-set protocols
- Up to 50 user-defined protocols

TECHNICAL SPECIFICATIONS OF THE BTL-CPMOTION™

Model	BTL-CPMOTION™ K EASY	BTL-CPMOTION™ K PRO	BTL-CPMOTION™ K ELITE
Part number	P093.001	P093.002	P093.003
User interface	START/STOP remote control	Remote control with touch screen	Remote control with touch screen
Advanced therapy settings		•	•
Preset protocols		•	•
User-defined protocols		•	•
Ankle module		○	•
Ankle extension/flexion			25°/40°
Knee extension/flexion		-10°/123°	
Max. movement speed		380°/min	
Dimensions		390 × 430 × 970 mm	
Weight	11 kg	12 kg	14 kg
Mains supply		100-240 V AC, 50-60 Hz	
Standard accessories	START/STOP remote control	Remote control with touch screen	Ankle module, remote control with touch screen

• Standard | ○ Optional